

General considerations of implant surgery

Körmöczi Kinga

Semmelweis University, Faculty of Dentistry



Department of Oral and Maxillofacial Surgery

**The planning of implant
rehabilitations happens
with the priority of the
prosthetic aspects, but in
close harmony with the
surgical ones.**



Treatment plan in implantology (backward planning)

- Protetical planning:
 - Type of toothless
 - Type of the prothesis
 - The protetical possibility for transfer the force optimally
- Sugical planning:
 - Place of implants
 - Number
 - Size, shape
 - Date of implantation


Model analysis before the surgery





The ideal mesio-distal position of implants



The ideal positioning of implants



Fabrication of X-ray, surgical template




Clinical aspects: Implant placement

- **Sterility**
- Speed of drilling
- Cooling
- Surgical technique



Clinical aspects: Implant placement

- Sterility
- **Speed of drilling**
- Cooling
- Surgical technique




Clinical aspects: Implant placement

- Sterility
- **Speed of drilling**
- Cooling
- Surgical technique

**Report of bone cells death
at temperature 40° C**

/Eriksson AR.,Albrektsson J.:
J.Prosth.Dent. 50:101,1983/



Surgical unit for implant placement:

- Adjustable speed and torque control
- Integrated cooling system
- Surgical hand-piece with cooling



Clinical aspects: Implant placement

- Sterility
- Speed of drilling
- **Cooling**
- Surgical technique




Outer, 5° C sterile
NaCl solution



Clinical aspects: Implant placement

- Sterility
- Speed of drilling
- Cooling
- **Surgical technique**

Intermittent drilling







Surgical steps of implant placement

- **Flap design**

- Marking drill
- Pilot drill
- Control measurements
- Widening drill
- Countersinking
- Thread cutting
- Implant placement

Various flap designs




**The flapless surgery can be
predictable and safe at ideal
anatomic conditions and proper
diagnostic evaluations.**

/Campelo L.D., Dominguez Camara J. R.,

Int. J. Oral Maxillofac. Impl. 2002; 17: 271-276/

Flapless implant placement /Straumann /



Surgical steps of implant placement

- Flap design
- **Marking drill 1200-1500 rpm**
 - Pilot drill
 - Control measurements
 - Widening drill
 - Countersinking
 - Thread cutting
 - Implant placement



Marking drill



Ø2.3mm




Ø3.1mm




Surgical steps of implant placement

- Flap design
- Marking drill
- **Pilot drill 1200-1500 rpm**
 - Control measurements
 - Widening drill
 - Countersinking
 - Thread cutting
 - Implant placement

Pilot drill /Ø 2.2 mm/



Ø2.2mm



Surgical steps of implant placement






- Flap design
- Marking drill
- Pilot drill
- **Control measurements**
 - Widening drill
 - Countersinking
 - Thread cutting
 - Implant placement

Controls

Depth gauge

Drill guide

Parallel gauge









Controls

Depth gauge

Drill guide

Parallel gauge





Controls

Depth gauge

Bur guide



Parallel gauge



Surgical steps of implant placement


- Flap design
- Marking drill
- Pilot drill
- Control measurements
- **Widening drill 600-800 rpm**
 - Countersinking
 - Thread cutting
 - Implant placement

Widening drill / \varnothing 2.8 mm/




\varnothing 2.8mm

Surgical bone condensation with Osteotom




The technique of bone condensation with Osteotom






Surgical steps of implant placement

- Flap design
- Marking drill
- Pilot drill
- Control measurements
- Widening drill
- **Countersinking** /case dependent/
 - Thread cutting 300 rpm
 - Implant placement

Geometry of the profile drill



Countersinking, with profile drill





2.8mm

Surgical steps of implant placement

- Flap design
 - Marking drill
 - Pilot drill
 - Control measurements
 - Widening drill
 - Countersinking
 - **Thread cutting** /case dependent/
 - Implant placement
- 25 rpm

Thread cutting





- Depends on the bone quality and the implant design







Surgical steps of implant placement

- Flap design
- Marking drill
- Pilot drill
- Control measurements
- Widening drill
- Countersinking
- Thread cutting
- **Implant placement**



Implant placement I.






Implant placement II.



Implant placement III.



Placement of cover- screw and suturing





**Healed
implant
before
second
surgery**