

APPLICATION FORM

SEMMEIWEIS UNIVERSITY FACS CORE FACILITY
BECKMAN COULTER CYTOFLEX SRT
1094 BUDAPEST, EOK, TÚZOLTÓ STREET. 37-47.

Flow cytometry analysis

Cell sorting

ORDERING INFORMATION

Name:

Email address:

Phone number:

Planned date:

Planned duration (hours):

Group Leader (person who should pay):

SAMPLE INFORMATION

Cell type:

Origin:

Description of risk (at least one week after viral infection):

Number of samples:

Estimated sample size (cell number for sorting):

Type of collection devices (5 ml, 15 ml tube or plate type):

List of fluorochrome/ fluorescent protein/ fluorescent dyes:

PROJECT DESCRIPTION

CYTOMETER SETTINGS

Temperature (5-37 °C):

Channels used:

405 nm Violet laser:	450/45
	525/40
	610/20
	660/10
	780/60
488 nm Blue laser:	525/40
	690/50
561 nm Yellow laser:	585/42
	610/20
	675/30
	710/50
	780/60
638 nm Red laser:	660/10
	712/25
	780/60