

Minimally invasive and micro-invasive dentistry. Digital dentistry.

Dr. Tamás Demeter

Semmelweis Universität Fakultät für Zahnheilkunde
Lehrstuhl für Propädeutik



SEMMELWEIS
UNIVERSITY 1769

Minimal intervention dentistry (MID)

- A **growing trend** in all fields of dentistry
- This term has been introduced by Davis and Makinson (1992)
- Instead of using the primary surgical and operative methods, this new concept is based on **four modalities**:
 - The early detection of caries and risk factors
 - Reduction of cariogenic microbes
 - Remineralisation of early lesions
 - Reparation and minimally invasive treatment of the evolved lesions
- The **main principle**: to keep all teeth and those functions through the life



Minimal intervention dentistry (MID)

This principle involves **more specialities** in dentistry

- Preventive dentistry
- Conservative dentistry
- Prosthodontics
- Orthodontics and Pedodontics
- Periodontology

Minimal intervention dentistry (MID) – in preventive dentistry


- **Early detection of** caries → Cariological diagnostics
- **Caries risk assesment**

- **Caries prophylaxis** with the propagation of natural remineralisation of enamel
 - Diethetic Control (change of food intake)
 - Plaque Control
 - Salivation Control
 - Instructions and motivation of the patient

- **Remineralisation** of the begunned lesions
- **Biofilm-modulation**
- Preventive closure of the caries predilection areas: **fissure sealing**



Caries diagnostics




- **Visual and tactile caries diagnostics**

- Probing – carefully, with blunt-ended probe!
- Limited diagnostic possibilities in the interproximal areas

- **X-Ray Diagnostics**

- Occlusal („Bitewing”) recordings
- Min. 60% loss of hard tissues!



Caries diagnostics

- Use of special devices:

- **Transillumination**

- Fiber-optics: FOTI (Fiber-optic Transillumination), DIFOTI
- Infrared: **DiagnoCam®**

- **Fluorescence based Systems**

- Infrared Laserfluorescence: **DiagnoDent, DiagnoDent Pen®**
- QLF (Quantitative Light-Induced Fluorescence) – Fluorescence based Intraoral camera: **VistaCam, VistaProof, SoproLife®**

- **Electric Impedance-based Systems**


- Electric Impedance-Spectroscopy: **CarieScan Pro®**

- **Photothermal Radiometry**

- A laser beam directed onto the tooth will be rebounded with heat production and luminescence → this will be perceived
- **Canary System®**



SOPROLIFE
Light Induced Fluorescence Evaluator



Remineralisation – Local (topical) use of Fluorides



- Fluoride varnishes
- Fluorid Iontophoresis
- Individual oral hygiene: toothpastes, mouthrinses



Biofilm Modulation

- **CPP-ACP Technology**

- Gc Tooth mousse ®
- Recaldent ®



Minimal intervention dentistry (MID) – in the conservative dentistry – use of special operative techniques

- Microinvasive Techniques
 - Cariesinfiltration
 - Ozone Therapy
- Minimalinvasive Techniques
 - Enamel microabrasion
 - Air abrasion, Hydroabrasion
 - Oscillation-based Instruments: Ultrasonic abrasive, Sonoabrasive Preparation methods
 - Laser Preparation
 - Rotary Instruments → Special cavity designs and Preparation techniques (Micropreparation; box-preparation; Tunnel-Preparation...)

Minimal intervention dentistry (MID) – in the conservative dentistry – use of special operative techniques

- **Microinvasive Techniques**

- Cariesinfiltration
- Ozone Therapy


- **Minimalinvasive Techniques**

- Enamel microabrasion
- Air abrasion, Hydroabrasion
- Oscillation-based Instruments: Ultrasonic abrasive, Sonoabrasive Preparation methods
- Laser Preparation
- Rotary Instruments → Special cavity designs and Preparation techniques (Micropreparation; box-preparation; Tunnel-Preparation...)

Caries infiltration



- Special microinvasive operative technique, which stops the enamel caries lesions and also the early stages of dentine caries, without irreversible removal of hard tissues
- **Principle of Caries infiltration:**
 - penetration of a **methacrylat-based, low viscosity material** (infiltrant) into the body of a carious lesion
 - After light-curing, the **infiltrant closes the micropores of the lesion**, and forms a barrier which prevents bacteria and substrates to diffuse into the lesion
 - the carious process will be **stopped**
- **DMG Icon ®**
- In our Phantom Head-Practice this therapy method will be explained and introduced



Caries infiltration - Steps

- tooth cleaning with dental floss and polishing paste
- absolute isolation with dental dam
- Inspection of the Lesion: is there a cavitation?
- Separation of the teeth with a special wedge (Icon; DMG, Hamburg)
- use of a special celluloid matrix with the etchant syringe (Icon Etch)
- Application of the 10 % hydrochloric acid (Icon Etch) for 2 minutes
- Rinsing off the acid with air-water syringe, then drying the surface
- Application of Ethanol (Icon Dry)
- thorough drying
- use of a new celluloid matrix with the infiltrant syringe (Icon Infiltrant)
- 1. Infiltration: 3min
- Removal of the surplus material with dental floss and air syringe
- Light polymerisation for 40 s
- 2. Infiltration: 1min.
- Removal of the surplus material with dental floss and air syringe
- Light polymerisation for 40 s
- Removal of the surplus material with scaler and finishing strips
- documentation!

Caries infiltration - Steps

Preparations

- Tooth cleaning with dental floss and polishing paste
- Absolute isolation with dental dam
- Inspection of the Lesion: is there a cavitation?
- Separation of the teeth with a special wedge (Icon; DMG, Hamburg)



Caries infiltration - Steps

Etching the surface

- Use of a special celluloid matrix with the etchant syringe (Icon Etch)
- Application of the 10 % hydrochloric acid (Icon Etch) for 2 minutes
- Rinsing off the acid with air-water syringe, then drying the surface
- Application of Ethanol (Icon Dry)
- Thorough drying



Caries infiltration - Steps

Application of the Infiltrant

- use of a new celluloid matrix with the infiltrant syringe (Icon Infiltrant)
- 1. Infiltration: 3min
- Removal of the surplus material with dental floss and air syringe
- Light polymerisation for 40 s
- 2. Infiltration: 1min.
- Removal of the surplus material with dental floss and air syringe
- Light polymerisation for 40 s
- Removal of the surplus material with scaler and finishing strips
- documentation!



Caries infiltration – in the esthetic region



Caries infiltration – in the esthetic region



Caries infiltration – in the esthetic region



Caries infiltration

Contraindications

- deep caries (ICDAS 4, 5 und 6 and/or radiological > D1)
- inactive carious Lesions
- Root caries (Dentin)
- Erosions


Restrictions

- carious Lesions in Pits and Fissures
- Molar-Incisor-Hypomineralisation syndrome (MIH), Amelogenesis imperfecta
- Fluorosis
- enamel development defects, traumatic enamel lesions

Ozone therap



- Non- or Microinvasive possibility to treat early fissure caries, without drilling and pain
- Endodontic use is also possible
- Mechanism of action: the tooth is treated with ozone gas for 20-60 seconds → bactericid, virucid, fungicid effect; The carious process will be stopped
- After treatment the treated tooth should be remineralized (fluoride varnish)
- The device is very expensive; therefore this treatment modality gained little popularity



Minimal intervention dentistry (MID) – in the conservative dentistry – use of special operative techniques

- Microinvasive Techniques

- Cariesinfiltration
- Ozone Therapy

- Minimalinvasive Techniques

- Enamel microabrasion
- Air abrasion, Hydroabrasion
- Oscillation-based Instruments: Ultrasonic abrasive, Sonoabrasive Preparation methods
- Laser Preparation
- Rotary Instruments → Special cavity designs and Preparation techniques (Micropreparation; box-preparation; Tunnel-Preparation...)

Enamel Microabrasion




- Removal of the surface stains of the enamel layer (until 0,2 mm depth)
- Clinically usable in cases of enamel hypocalcification, hypoplasty, fluorosis or demineralisation
- Principle: decalcification of the enamel (with 6% hydrochloric acid), then removal of a 20-25 micrometer thick layer from the surface (with SiC Abrasive Paste)
- This treatment can also be completed by adhesive filling therapy in the operative field
- Opalustre ®



Air abrasion; Hydroabrasion




- Minimalinvasive cavity preparation technique – without drilling
- Mechanism of action: blowing of 27-50 micron Aluminiumoxid-Particles with air pressure onto the tooth surface
- Use: preparation for Fissure sealing; Minimalinvasive treatment of Fissure caries; surface roughening for adhesive procedures
- PrepMaster®, EtchMaster®, Kavo Rondoflex® (Hidro-abrasion)





Ultrasonic / Sonoabrasive Preparation

- ultrasonic or sonic (25 kHz/6,5 kHz), pneumatic preparation instruments
- Various working heads for each cavity design
 - endodontics, Inlay-Preparation, Approximal cavities
- Controlled tooth structure removal, with maximum protection of the neighbouring teeth
- **Kavo SonicFlex®**



Laser Preparation

- Possibility of selective tooth preparation and Caries removal - without Vibration und Sensibility
- After preparation the surface will be rough, with opened dentinal tubules → good for adhesive techniques
- Mechanism of action: The carious tissues with high water content will be removed with ablation
- Erbium Laser (with Infrared Wave length) is the best for cavity preparation and caries removal: Er-YAG; Er-YSGG
- Waterlase ®



Use of Rotary Instruments new way - "Micro-dentistry"



- The original rotary instruments can also be used with the principles of minimally invasive dentistry
- Use of special burs, which guarantee a controlled and minimally invasive reduction of tooth structure: **MicroPrep Kit (KOMET®)**
- Use of the Operative microscope and various loupes in the modern dentistry:
 - Cavity preparation
 - Endodontics
 - Minimally invasive Prosthodontics: Preparation and cementing




Special cavity designs

- Approximal Box Preparation
- Tunnel Preparation




Digital Dentistry

The basics & possibilities



CAD/CAM Systems

- Digitising device (Intraoral scanner)
- Data processing software (CAD - **C**omputer-**A**ided **D**esign)
- Manufacturing technology (CAM - **C**omputer-**A**ided **M**anufacturing)
 - Subtractive techniques (Dental milling machine)
 - Additive (building) techniques (3D printer)
- To complement or replace traditional processes.




From the teaching material of Dr. Zsófia Vincze

Digitising device

- **Optical scanner**

- "Triangulation" - compiles three-dimensional data series based on the angle between the light source and the receiver



3D scanned data file in STL format

- **Mechanical scanner**

- The pattern is traced line by line using a round head part




From the teaching material of Dr. Zsófia Vincze

Intraoral scanner



From the teaching material of Dr. Zsófia Vincze

Planning software



- Design of restorations (crown, bridge, primary telescope, etc...)
- STL file
- Continuous improvements, updates available

From the teaching material of Dr. Zsófia Vincze




Dental milling machine

From the teaching material of Dr. Zsófia Vincze

Dental milling machine

- CNC (Computer Numerical Control) milling machine
 - Programmable
 - Contains a computer itself – controls the movement of the milling tool
- CAD designed restoration (e.g. bridge framework)
 - Needs to be complemented with holding pins
- STL – initial data for CAM program
 - Planning the spatial movement of the milling tool
 - Data sent to CNC milling machine control computer



From the teaching material of Dr. Zsófia Vincze


3D Printing



From the teaching material of Dr. Zsófia Vincze


3D Printing

- Additive technology - adding model material layer by layer
- CAD designed model is broken down into slices - layer by layer the 3D object is built
- Problems encountered: Material, strength, stiffness, unsupported parts, size, accuracy...



From the teaching material of Dr. Zsófia Vincze

3D Printing



From the teaching material of Dr. Zsófia Vincze

Making a dental bridge

Traditional vs digital technology



From the teaching material of Dr. Zsófia Vincze

Traditional



Digital




From the teaching material of Dr. Zsófia Vincze

Traditional



3D printed

Digital





From the teaching material of Dr. Zsófia Vincze



Traditional




Digital



From the teaching material of Dr. Zsófia Vincze

Traditionell



From the teaching material of Dr. Zsófia Vincze

Back Preview Next

Anatomy design

Temporary on prepared model 23

Smile Composer

Smile Library 2

Sculpt

Multiple teeth tools

Single tooth tools

Smart tools

Connection to margin line

Wax knife

2.16 mm

0.80

1 2 3 4 5 6 7

Sculpt protected surfaces

Lift anatomy off

Digital



Additional scans

Load Align Save

Save model Keep additional information

Please select a model to save. Hold Ctrl to add multiple models for saving in one file.

Order: Temporary on prepared model 13, 12, 11, 21, 22, 23
 Clinic: Rafael Galakberow

From the teaching material of Dr. Zsófia Vincze



Dental milling machine



3D Printing



From the teaching material of Dr. Zsófia Vincze



From the teaching material of Dr. Zsófia Vincze

Thank you for your attention



SEMMELWEIS
EGYETEM 1769